



# **Product specifications**

Name	Anti-h SAA 2211 SPTN-5				
Specificity	Antibody recognizes human serum amyloid A (SAA)				
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.				
Product code	100808				
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative				
Shelf life and storage	Unspecified, storage at 2–8 °C				
Subclass	IgG <sub>1</sub>				
Analyte description	Serum amyloid A, a family of apolipoproteins, is associated with high density lipoprotein during inflammatory states. The level of serum amyloid A (SAA) in the blood increases dramatically in response to tissue injury and inflammation. SAA also acts as a cytokine, influencing cell adhesion, migration, proliferation and aggregation. SAAs are implicated in several chronic inflammatory diseases, such as amyloidosis, atherosclerosis, and rheumatoid arthritis.				

# Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80–120 % compared to the reference sample in an FIA test
IEF Profile	6.3–6.7
Purity	≥ 95 %

## **Kinetic parameters**

Association rate constant	4.9 x 10 <sup>5</sup> 1/Ms
Dissociation rate constant	2.5 x 10 <sup>-₅</sup> 1/s
Affinity constant	K <sub>A</sub> = 2.0 x 10 <sup>10</sup> 1/M K <sub>D</sub> = 4.0 x 10 <sup>-11</sup> M (= 0.04 nM)
Determination method	BLI (Octet RED96e)
Determination antigen	Recombinant SAA, Medix Biochemica (Cat 610070)

### **Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532 Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.





### **Cross-reactivities**

Epitope

Recognizes human SAA1, human SAA2, canine and equine SAA. Does not recognize human SAA4, feline or bovine SAA.

#### Conformational

FIA, LF, IT

### Pair recommendations

		DETECTION								
		2201/ RC2201	2203	2205	2208	2209	2211	2212	2213	
CAPTURE	2201/ RC2201	-	+	+	-	-	+	-	-	
	2203	+	-	+	-	+	+	-	+	
	2205	-	-	-	-	-	+	+	-	
	2208	-	-	-	-	-	-	+	-	
	2209	-	+	-	-	-	+	-	-	
	2211	+	-	-	-	+	-	-	-	
	2212	-	-	+	+	-	-	-	-	
	2213	-	-	-	-	-	-	-	-	

Following pairs are especially recommended for the below mentioned assays: FIA: 2205 (capture) – 2212 (detection), 2208 – 2212 LF: 2201 (membrane) – 2203 (particles), 2203 – 2201 IT: 2201 – 2211, 2201 – 2203, 2209 – 2211

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Native human SAA, EastCoast Bio (Cat LA449), recombinant human SAA, Medix Biochemica

Platforms tested

Antigens tested

Product stability

 TEMPERATURE, TIME
 RESULT

 -70 °C, 21 days
 OK

 -20 °C, 21 days
 OK

 +4 °C, 21 days
 OK

 +35 °C, 21 days
 OK

 +45 °C, 7 days
 OK

(Cat 610070) and Lee Biosolutions (Cat 527-18).

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

#### Miscellaneous

References

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532 Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.